



SAFETYBOX® AW 120 NEW

Data Sheet

SAFETY CABINET FOR THE STORAGE OF 160 LITRES OF CHEMICAL, ACID AND BASE PRODUCTS

- Completely produced in pickled steel, 10/10 mm thickness, painted with epoxide acid-proof dusts and passage through thermic tunnel at 200° C. Thickness of painting from a minimum of 60µm (micron) to a maximum of 250µm (micron) as per UNI EN ISO 9001:2008 procedure.
 - **Certified palletized base:** new design of the base for easier transportation of the cabinet. **With removable cover.**
 - **External monolithic construction** completely sealed. New building method with double cases, internal and external.
 - **3 shelves adjustable in height** made of electrogalvanized steel painted with epoxide acid-proof dusts. Tray-shaped shelves to collect any liquid spillage of breakages.
Removable and adjustable shelves with anti-falling security lock by rack.
 Collection capacity of the shelf approx. 15 litres. Loading capacity kg. 100.
 - Watertight bottom basin in electrogalvanized steel painted with epoxide anti-acid dusts. Capacity approx. 44 litres.
 - Doors opening up to 110° to allow an easy extraction of the shelves without tilting them.
 - **Window doors with unbreakable glass.**
 - Security lock with key, multi point closure and extractable cylinder.
 - Acid-proof hidden hinges.
 - “Danger” warning stickers on the doors (as per European Norms DIN 4844 and ISO 3864).
 - Grounding connection to discharge build-up of static electricity.
 - Levelling feet for balance.
 - **Lifetime warranty “RC PRODOTTI”:** it covers any injury to persons and objects.
 - Usage and maintenance manual.
 - **Packaging with “Tipping over system®” device.**
-
- External dimensions: mm 1200 x 500 x 1690 h
 - Internal dimensions: mm 1120 x 464 x 1522 h
 - Weight: Kg 140
 - Internal volume: approx. 824 litres

ACCESSORIES

- Possibility to provide the cabinet with aspirating and filtering kit (KIT 01)
- Additional shelves